



Healopad Hardware Specifications:

Dimensions: 11 X 6X 3 cm

Weight: 120g

Body Material: Plastic

Cable Length: Can be up to 1 meter long

Battery Charger: 12V DC Adaptor, One time full charge works for 7 days minimum.

Not Waterproof.

Please keep the device clean and dry.

The BOX Unit Appearance Guide

- CONN BT: It shows connectivity. If not connected will blink per second. If connected will blink once, each 3 seconds.
- BAT STAT: When the battery needs charging, it will be in Red.
- CHARG CONN: A 12V of DC adaptor can charge the unit.
- SN: Shows your unit serial number.

Healopad Smart Wound Bandage Box User Manual:

- Install the App on your mobile.
- Login
- Turn the power on.
- Click Scan in the App and choose BLE in the shown list.
- On the left upside of the page, “Connected” will appear if everything is ok and data transmission will be started.